QuikRead go [®] HbA1c Control Set	Cat. no. 154520	English
<div>Intended use</div> <p>QuikRead go[®] HbA1c Control Set includes controls with quantitative assigned values for monitoring the performance of quantitative QuikRead go[®] HbA1c assay used with the automated QuikRead go[®] instrument. The control set is intended to be used by trained healthcare professionals in clinical laboratories and near-patient testing (NPT) settings. For <i>in vitro</i> diagnostic use.</p> <div>Principles of the procedure</div> <p>Use the controls in the same way as unknown blood specimen would be used, following the instructions for the assay procedure.</p> <div>Reagents</div> <p>The QuikRead go HbA1c Control Set features two HbA1c controls with different HbA1c concentrations: QuikRead go HbA1c Control and QuikRead go HbA1c Control High. The control is human-based material. It should be handled in the same way as patient samples. The control is stable for 5 minutes in the cuvette with a closed cap.</p> <p>Content of the kit:</p> <ul style="list-style-type: none">- 1 x 0,5 ml QuikRead go HbA1c Control - 1 x 0,5 ml QuikRead go HbA1c Control High - Instructions for use <div>Preparation</div> <p>The control is ready-to-use. Mix well. Do not shake or turn upside down to avoid contact with the rubber stopper. Keep the rubber stopper closed during the storage.</p> <div>Storage</div> <p>Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 50 days at 2...8°C. Do not use any vial that has visible evidence of microbial growth.</p> <div>Expected values</div> <p>The individual results for the Control and Control High should be within the range indicated on the label of the control package for each control respectively. The target values and acceptable range for the controls are reported in mmol/mol HbA1c (aligned to the IFCC</p>		

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Reference Procedure^{1,2}) and in % HbA1c (aligned to the DCCT study assay, also known as NGSP-HbA1c^{3,4}). Each laboratory should establish its own mean and precision data. The values are traceable to IFCC Reference Measurement Method¹.

Warnings and precautions

- Reagent contains <0,1% tetrachloroethylene, which is not considered a harmful concentration. Tetrachloroethylene may cause eye irritation and slight skin irritation. In the case of contact with eyes, rinse immediately with plenty of water for 15 minutes and seek medical attention. If contact with skin, wash off with soap and plenty of water. Consult a physician. Thermal decomposition of tetrachloroethylene can lead to release of irritating/toxic gases and vapors. Tetrachloroethylene reacts with strong oxidizers and chemically active metals such as barium, lithium, and beryllium.
- Avoid using the last drops of control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new vial.
- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- The control material is sensitive to temperature. Avoid unnecessary storage outside 2...8°C.
- The control contains material of human origin. All blood donors have been tested and found negative by FDA approved methods for HbsAG, HIV-1, HIV-2 and HCV. As a precaution, handle the control as if capable of transmitting infectious agents.

- Disposal**
- Disposal of contents according to national and local law.
 - Materials of the components
 - Paper: Instructions for use
 - Cardboard: Outer box
 - Glass: Control vials
 - Plastic: Vial caps
 - Rubber: Stopper
 - When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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