

Convenient testing of several persons simultaneously

Pulsar Ecotest COVID-19 Antigen
Saliva Test



Fast COVID-19 result



AIDIAN

Pulsar Ecotest COVID-19 Antigen Saliva Test

The rapid and convenient solution for detection of COVID-19 infection in point of care.

Pulsar Ecotest COVID-19 Antigen Saliva Test's innovative design enables safe handling and convenient use and makes it especially **handy for testing several people simultaneously**. One single packed test contains all materials necessary for testing making the distribution of tests to separate testing sites easy. The test is also simple to take along to the actual site of testing, when needed.

Pulsar Ecotest COVID-19 Antigen Saliva Test detects the SARS-CoV-2 nucleocapsid protein from saliva samples.


Aidian Oy

aidian@aidian.eu / +358 10 309 3000

For orders: orders@aidian.eu

Pulsar Ecotest COVID-19 Antigen Saliva Test, 20 tests

Product no COV-P53	Ordering no 154860
-----------------------	-----------------------



Fast and safe

- Only a few steps to the test result
- Hands-on time just seconds
- Closed system minimizes contamination risk
- Enables fast testing of several people simultaneously

Easy and convenient

- Pleasant sampling
- Prevents user errors
- Single packed tests allow easy distribution
- Storage at room temperature

Reliable results

- Detects SARS-CoV-2 N-proteins (SARS-CoV-2 variants are found in S-protein)
- Compared to PCR Relative sensitivity 90%, specificity 97%
For Ct <25, Sensitivity 98.3%, Specificity 100%

For *in vitro* diagnostics.

The test result is used for early detection of a possible infection by medical professionals. A confirmation test is necessary.